
Quadra Electrophoresis Systems

4-Place Mini-Vertical Slab Gel /Blotting Systems

QNC-700, with cooling

QNX-700

INSTRUCTION MANUAL

T A B L E O F C O N T E N T S

Important User Information	3-4	
Section 1: General Information		
1.1 Introduction and Features	5	
1.2 Specifications	5	
1.3 Comb Specifications for Hand Cast Gels	6	
1.4 Safety	6	
Section 2: Description of Parts		
2.1 Unpacking and Components	7	
Section 3: Instructions for Slab Gel Electrophoresis Using Pre-Cast Gels		
3.1 Preparing Electrophoresis Unit Using Pre-Cast Gels	8-9	
3.2 Running the Gel	10	
3.3 Removing the Gel	11	
Section 4: Alternate Protocol for Slab Gel Electrophoresis Using Hand-Cast Gels		
4.1 Gel Wrap™ Gasket Casting Method	12-13	
4.2 Preparing Electrophoresis Unit Using Hand-Cast Gels	14	
Section 5: Instructions for Western Blotting		
5.1 Preparing Unit for Electro-Blotting	15-16	
5.2 Electro-Blotting Procedure	16	
5.2 Removing the Blot	16	
Section 6: Running Conditions		
6.1 Recommended Power for Slab Gels	17	
6.2 Recommended Power for Electro-Blotting	18	
6.3 References	18	
6.4 Recommended ClearPAGE™ Buffer Formulation	19-20	
Section 7: Maintenance of Equipment		
7.1 Care and Handling	21	
7.2 Maintenance	21	
Section 8: Ordering Information		
8.1 Ordering Information for Quadra Electrophoresis System	22	
8.2 Ordering Information for Related Products	23	
Section 9: Contact Information		24

IMPORTANT USER INFORMATION

This Instruction Manual will explain how to use this product safely and effectively. Please read and carefully follow the instruction manual in its entirety.

The triangle/exclamation mark symbol alerts the user of the product to important operational, maintenance, and/or warranty requirements.

The triangle/lightning bolt symbol alerts the user of the product to potentially hazardous electrical exposure.

Failure to adhere to the instructions could result in personal and/or laboratory hazards, as well as invalidate any warranty. Always turn off the DC power source prior to disconnecting power cords from the product. Disconnect power cords from the power source first and then from the product. For maximum safety, always operate this system in an isolated, low traffic area, not accessible to unauthorized personnel. Never operate damaged or leaking equipment.

WARRANTY AND LIABILITY

This product was produced utilizing the highest practical standards of materials, workmanship, and design. C.B.S. Scientific warrants that the product has been tested and will meet or exceed published specifications. This warranty is valid only if the product has been operated and maintained according to the instructions provided.

C.B.S. Scientific warrants this product to be free from defects in materials and workmanship under normal service for one year from date of shipment. If the product proves defective during this period, C.B.S. Scientific will repair or replace it at our option, free of charge, if returned to us postage prepaid. This warranty does not cover: damage in transit, damage caused by carelessness, misuse or neglect, normal wear through frequent use, damage caused by solvent corrosion, damage caused by improper handling or user alteration, nor unsatisfactory performance as a result of conditions beyond our control. C.B.S. Scientific shall in no event be liable for incidental nor consequential damages, including without limitation, lost profits, loss of income, loss of business opportunities, loss of use and other related damages, however caused, nor any damage arising from the incorrect use of the product.

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General Information

1.1 Introduction

C.B.S. Scientific offers the Quadra Electrophoresis Systems for performing SDS-Page, acrylamide-nucleic acid separations and electro-blotting. These quadruple capacity units provide the capability of running or blotting four gels simultaneously. The Quadra comes in two models, with (QNC-700) or without (QNX-700) temperature controlled buffer conditions. Each unit includes 2 cores, 4 blotting cassettes with sponge pads.

These systems are designed primarily to be used with a variety of precast gels. (See Table 1). Hand-cast gels can also be run in this unit with the use of accessory combs, plates and spacers.

Table 1: Features of Vertical Mini-Gel Systems

Model #	Capacity	Plate Dimensions (w x h)	Compatible Precast Gels	Cooling
QNC-700 QNX-700	Quadruple	10cm x 10cm 10cm x 9cm* 10cm x 8cm *requires adaptor shim plate	ClearPAGE Expedeon Invitrogen PAGER NuPAGE SERVAgel Lonza	QNC-700 has double-sided cooling labyrinth for use with external cooler

1.2 Specifications

Constructions:

Buffer chamber, safety cover, core

Electrodes

Power cords

Combs

Glass plates

Spacers

Clamps

Safety Certification

Acrylic

Pure Platinum wire .009" diameter

FR Silicone rated 7500VDC, 200mA, 65°C

Teflon

Soda-lime float glass & Borosilicate

PVC

Polypropylene, stainless steel

EN61010-1-1993 (IEC1010-1)

Table 2: Specifications

Model #	QNC-700, QNX-700
Shipping Weight	10lbs
Overall Size	40cm(l) x 30cm(w) x 25cm(h) 16"(l) x 12"(w) x 10"(h)
Recommended buffer volume: a) Separation phase: <ul style="list-style-type: none"> • 2 cathode reservoirs • 1 anode lower reservoir • Total Buffer min./max. requirement 	200mls/each for a total of 400mls 1100mls 0.95 - 1.5 L
b) Electro-blotting phase	1.7 L total
Distance between electrodes	9cm
Voltage limit	250V

1.3 Comb Specifications for Hand Cast Gels

COMB OPTIONS, material: PTFE, tooth depth: 9.5mm, overall length: 7.9cm				
Cat. #	# of Teeth	Thickness of Teeth (mm)	Width of Teeth (mm)	Recommended Max. Sample Vol./Well (µl)*
MVX-0701	1	0.75	66.6	356
MVX-0710	10	0.75	5.23	28
MVX-0714	14	0.75	3.3	19
MVX-0716	16	0.75	2.92	16
MVX-1001	1	1.0	66.6	474
MVX-1010	10	1.0	5.23	37
MVX-1014	14	1.0	3.3	24
MVX-1016	16	1.0	2.92	21
MVX-1501	1	1.5	66.6	712
MVX-1510	10	1.5	5.23	56
MVX-1514	14	1.5	3.3	35
MVX-1516	16	1.5	2.92	31

*Maximum loading volume is calculated by taking 75% of the total tooth volume.

1.4 Safety

Power to the Quadra Electrophoresis System is to be supplied by an external DC voltage power supply that must be ground isolated so that the DC voltage output floats with respect to ground. For any power supply used, the maximum specified operating parameters for the units are:

Maximum Limits

300 VDC

30 watts power

190mA current

60°C ambient temperature

Current to the unit, provided from the external power supply, must enter the unit through the lid assembly, providing a safety interlock to the user. Current to the unit is broken when the lid is removed. **Do not attempt to use the unit without the safety lid, and always turn the power supply off before removing the lid, or when working with the unit in any way. Follow safety precautions specified by the power supply manufacturer.**

SECTION 2


Description of Parts

2.1 Unpacking and Components

Please verify that your unit comes complete with the following components:

Quadra Electrophoresis System:

- Lower Reservoir (QNC-700 includes 2 cooling labyrinths and 2 tubing adapters)
- Safety cover
- Power leads (must be attached before using. See section 3.1.1)
- 2 Cores
- 4 blotting cassettes with foam pads (not shown)
- Adaptor plate for running single gels (not shown)

SECTION 3

Instructions for Slab Gel Electrophoresis Using Pre-Cast Gels


3.1 Preparing the Electrophoresis Unit

1. Power leads are shipped separately and **must be attached before using your Quadra**. Thread the power lead ends into the lid receptacles and rotate in clockwise direction to hand tighten.


2. Place unit in authorized work area. Remove safety cover from the assembled unit by simultaneously pressing down on white push pins while lifting up on blue safety cover as shown in figure 1.


Do not remove safety cover by pulling up on leads!


3. Remove white core from lower reservoir by grasping core with one hand and lifting directly up as shown in figure 2. Remove second core in same manner.


4. Open doors on the core assemblies by pulling up on the white latches, as shown in figure 3.


5. Slide pre-cast gel cassette or plate set(s) into the core assemblies with the notched plate facing in towards the upper buffer reservoir as shown in figure 4. If using a pre-cast gel stored at 4°C, allow to warm to room temperature. If pouring your own gels please see Alternate Protocol (Section 4) on pages 12-14 describing gel casting using Gel Wrap™.

6. If running one gel, slide white plastic adaptor plate into the side without the gel. If running a Lonza 10cm x 9cm gel cassette place clear shim as shown in figure 5.

7. Close doors and relatch by pressing down on the white latches so that the assembly looks like that shown in figure 6.

3.2 Running the Gel

1. Place stirring bar into bottom of reservoir in stirring corral.





2. Place core assemblies into lower reservoir. The anode (red) and cathode (black) electrodes are color-coded on both the core/cassette assembly and lower reservoir. See figure 7. Ensure the red dot on the cassette assembly is on the same side as the red receptacle on the lower reservoir. Fill both core upper reservoir with freshly prepared buffer (~ 400mls).

If any buffer is spilled into banana jack receptacles (outlined in yellow boxes in figure 7) in lower reservoir, dry completely using compressed air! Failure to do this will result in accelerated banana jack corrosion. Quadra should look like figure 8.

3. Pour only enough freshly prepared buffer into lower chamber so that the final buffer level is just below bottom of sample wells. Using a pipette or syringe, thoroughly flush out the wells in the glass plate sandwich with buffer. Load samples. If outer lanes do not contain sample, it is recommended that you run standards and/or fill outer lanes with loading buffer to reduce smiling and wrap-around effects.

4. Attach safety cover and turn on magnetic stirrer. If you have the cooled unit, cat. #QNC-700, precool buffers and attach external cooler (set to the desired temperature) to the tubing adaptors so that the inlet is on the bottom and the outlet is on top as shown in figure 9. The Quadra is now ready to attach to the power supply.

5. Connect the leads to the power supply, matching the color-coded red to red and black to black as shown in figure 10. **See Section 6.1 for recommended power conditions.** Begin separation by electrophoresis.


3.3 Removing the Gel

1. **Turn the power supply off and disconnect the leads from the power supply.** Remove the safety cover from the unit, by placing thumbs on white posts next to red & black connectors, then pushing down while pulling up with fingers under lid as shown in figure 9. **Do not remove safety cover by pulling up on leads!**




2. Pull up on gel door latches, and open gel door. Remove gel sandwich from cassette assemblies. Stain and fix according to your preferred method.

4.1 Alternate Protocol: Using Gel Wrap™ Gasket Casting Method (Not using Pre-cast gels)

1. Place all components in an authorized work area. You will need: glass plate set, Gel Wrap™ Gasket, spacer set, comb, 3 GPC-0002 clamps, and polyacrylamide solution. Prepare and clean glass plates by hand washing both plates with a high quality lab detergent followed by a complete rinsing with dH₂O. Air-dry or use a lint-free tissue. Spray/wipe the chosen inner surfaces of the plate set with 95% ethanol and dry with lint-free tissue.
2. Start gel casting procedure by holding the 3mm thick, notched back plate with the rounded bottom corners and applying the gasket around one side of the glass plate. Note: one side of the “U” shaped gasket is flat, and the other side has tubing that will act as a seal around the spacers.
3. When applying the gasket over the rounded corners of the notched glass plate, make sure the cuts on the gasket align with the rounded corners of the glass plate. Once the gasket is pushed over the bottom edge and corners, work it down the remaining side.
4. Place the gasketed plate on the lab bench with the tubing side up, and extend the bottom of the plate over the edge of the bench, approximately $\frac{3}{4}$ of an inch. Place the spacers along side the inside edges of the gasket. Be sure the rounded corner end of each spacer is facing the outside bottom of the plate, following the radius of the glass.
5. Place the thinner unnotched back plate on top of the bottom assembly, starting from the bottom edge and gently easing the plate down. Verify the gasket is smooth around the edges and then clamp along the bottom.




6. Lift the assembly and stand it on the base of the clamp. For leveling, push glass plate assembly down until it stops against clamp body. Clamp the sides of the assembly with additional casting clamps on either side. As each clamp is attached, be sure the gasket is aligned between the plates forming a seal.

7. Apply PAGE solution to gel plate sandwich using a syringe or pipette. If using a stacking gel, pour desired height of running gel, then overlay a small amount of dH_2O or 0.1% SDS solution to top of gel.

8. After polymerization, rinse with buffer, add stacking gel solution and insert comb. For regular, unit percentage gels, add polyacrylamide solution to correct height, and insert comb. Allow gel to set, usually 20 minutes. Extra gel solution in pipette or syringe can be monitored to test polymerization of gel mix.

9. Disassembly (see right and below right). Hold the clamped plate assembly with one hand. Remove the gasket by starting at one of the top ends and pulling up and out on the gasket until it releases from the plate, up to the bottom of each of the white clamps. When each clamp is reached DO NOT remove it, instead feed the gasket down through the clamp body and repeat pulling up and out. Continue feeding until the gasket is fully detached. Once gasket is removed, detach clamps. If gel is not to be used immediately, wrap entire plate sandwich with plastic wrap tightly to seal and store at 4°C for up to a month.


SECTION 4

Alternate Protocol for Slab Gel Electrophoresis. (Not using Pre-cast gels)


4.2 Preparing the Electrophoresis Unit when Using Gel Wrap™ Gasket Casting Method



1. Place unit in authorized work area. Remove safety cover from the assembled unit by simultaneously pressing down on white push pins while lifting up on blue safety cover as shown in figure 1. **Do not remove safety cover by pulling up on leads!**

2. Remove white core from lower reservoir by grasping core with one hand and lifting directly up as shown in figure 2. Remove second core in same manner.
3. Open doors on the core assemblies by pulling up on the white latches, as shown in figure 3.
4. Slide glass plate sandwiches into the core assemblies with the notched plate facing in towards the upper buffer reservoir as shown in figure 4.
5. If running one gel, slide white plastic adaptor plate into the side without the gel.

6. Close doors and relatch by pressing down on the white latches so that the assembly looks like that shown in figure 5.
9. For rest of protocol please refer to instructions on pages 10-11 (Sections 3.2 and 3.3) for running the gel and removing it after electrophoresis.

SECTION 5

Instructions for Western Blotting


5.1 Preparing the Unit for Blotting

1. Remove safety cover from the assembled unit by simultaneously pressing down on white push pins while lifting up on blue safety cover. **Do not remove safety cover by pulling up on leads!** Remove white core from lower reservoir by grasping core with one hand and lifting directly up. Open doors on the core assembly by pulling up on the white latches, as shown in figure 1.


2. Open blotting cassette as shown in figures 2-3 and lay it flat on the bench.

3. Assemble blotting stack as shown in figure 4. With cassette wide open assemble components on black side in the following order: foam pad, gel*, buffer saturated transfer membrane, then buffer saturated blotting paper. Smooth with gloved finger or roll with glass rod to be sure no bubbles exist between the gel and the transfer membrane.


***Note:** to prepare gel for blotting, trim off wells and any excess acrylamide at the bottom, and invert 180° so that the large molecular weight proteins are at the bottom of the cassette. This puts them in contact with a stronger field strength and allows the blotting transfer to take place more efficiently.

4. Insert blotting cassettes into core making sure that red side faces outward. See diagram 5.

5. Close doors and re-latch by pressing down on the white latches so that assembly looks like that shown in figure 6. If running one blot, slide white reservoir conversion plate into the side without the blotting cassette.



5.2 Electro- Blotting Procedure

1. Place stirring bar in bottom corral of lower reservoir. Place core/blotting cassette assemblies into lower reservoir. The anode (red) and cathode (black) electrodes are color-coded on both the core/cassette assembly and lower reservoir. Ensure the red dot on the cassette assemblies are on the same side as the red receptacle on the lower reservoir. Attach safety cover and turn on magnetic stirrer. If you have the cooled unit, cat. #QNC-700, precool buffers and attach external cooler (set to the desired temperature) to the tubing adaptors so that the inlet is on the bottom and the outlet is on top (see fig. 9 on page 10).

2. Pour 1.7 liters of freshly prepared, chilled (4°) buffer into lower buffer reservoir. Buffer will percolate into central core.

3. Attach safety cover. The unit should look as shown in figure 7 and is ready for power.


4. Connect the leads to the power supply, matching the color-coded red to red and black to black as shown in figure 8. **See Section 6.2 for recommended power conditions.** Begin transfer by electrophoresis.

5.3 Removing the Blot

1. **Turn the power supply off and disconnect the leads from the power supply.** Remove the safety cover from the unit, by placing thumbs on white posts next to red & black connectors, then pushing down while pulling up with fingers under lid. **Do not remove safety cover by pulling up on leads!**

2. Blotting cassettes can be removed by leaving the core in place and opening the top latches of the core, opening the doors and lifting the cassettes out. Unlatch the blotting cassettes and remove blot from blotting sandwich.

SECTION 6

Running Conditions

6.1 Recommended Power for Slab Gels:

Precise electrophoresis conditions will vary according to the number and type of gels used, buffer conditions employed, power input, and the general goal of the experiment. Refer to reference section 6.3 for in depth discussions on practical and theoretical approaches to protein gel electrophoresis.

6.1.1 ClearPAGE Gels

Run Voltage	Starting Current	Ending Current	Approx. Run Time
180VDC	90mA/gel	40mA/gel	30-75 minutes

6.1.2 General Recommendations

If running only one gel, keep the volts the same but reduce the mA's by half. Keep in mind that as the thickness of gel increases, the mA's increase proportionally.

At constant voltage, the proteins will migrate at a constant rate during electrophoresis with adequate heating appropriate for denaturing gels. Increasing the voltage/ mA (for a single gel thickness and percentage) will speed mobility but increase the risk of overheating.

If using the cooling capability of the QNC-700 the power input and the migration rate can be increased. The joule heating generated by the higher power is offset by the cooling effect of the buffer between the gels.

6.1.3 Tris-Glycine Gels

For SDS-PAGE Tris-Glycine (Laemmli) buffer systems with **two** 1.0mm thick gels at room temperature use the following conditions at constant voltage:

80VDC until samples have fully entered stacking gel
120VDC @ 60mA-90mA/gel (depending on gel type) thereafter until dye is near bottom of gel

6.2 Electro-Blotting

As a general recommendation, equilibrate gels (after running) with the diluted transfer buffer for 5 to 10 minutes before transfer.

Blotting Buffer	
ClearPAGE™ Transfer Buffer 10X/20X cat. # FB82500 (see Section 6.4)	100ml (1:10 dilution)
Methanol	200ml
Ultrapure Water	720ml

Typical Blotting conditions for Quadra	
Power Supply Setting	200V constant
Blot time	1.5 - 2.0 hours with stirring
Expected current	180mA / 1 gel 220mA / 2 gels 300mA / 4 gels

6.3 References

1. Hames, B.D. (ed.) (1998). *Gel Electrophoresis of Proteins. A Practical Approach*. 3rd edn. Oxford University Press, Oxford. Ch. 1,3.
2. Sambrook, J., Fritsch, Russell, D. (2001). *Molecular Cloning. A Laboratory Manual*. 3rd edn. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, New York. A8.40-A8.55
3. Ausubel, F.M., Brent, R., Kingston, R.E., Moore, D.D., Seidman, J.G., Smith, J.A., Struhl, K. (ed) (1993). *Current Protocols in Molecular Biology*. Vol. 2, Greene Publishing Associates, Inc. and John Wiley & Sons, Inc., Ch. 10.

6.4 Recommended ClearPAGE™ Buffer Formulations

As an alternative to the ClearPAGE buffers available for purchase, these formulations may be used to prepare buffers yourself. Use high-quality, low-conductance ingredients. Do NOT use acid or base to adjust the pH!

Standard SDS Running Buffer, 20X for Reduced Samples (Cat. # FB60500)

* pH should be between 8.4 and 8.5 at 25° C.

Ingredient	MW	Molarity	Qty/Liter
Tricine (free acid)	179.17	0.8M	143.4 g
Tris (free base)	121.14	1.2M	145.2 g
SDS (2%)	288.38	-	20.0 g
Sodium Meta-bisulfite	104.06	50mM	5.0 g
Ultrapure water (fill to)	-	-	1000ml

* For non-reduced samples (especially antibodies), omit the Sodium Meta-bisulfite

Standard SDS Running Buffer, 20X for Non-Reduced Samples, (Cat. # FB50500)

* pH should be between 8.4 and 8.5 at 25° C.

Ingredient	MW	Molarity	Qty/Liter
Tricine (free acid)	179.17	0.8M	143.4 g
Tris (free base)	121.14	1.2M	145.2 g
SDS (2%)	288.38	-	20.0 g
Ultrapure water (fill to)	-	-	1000ml

* For non-reduced samples (especially antibodies), omit the Sodium Meta-bisulfite

Turbo SDS Running Buffer, 20X, (Cat. # FB14500)

* pH should be between 8.3 and 8.4 at 25° C.

Ingredient	MW	Molarity	Qty/Liter
MPS (free acid)	209.26	0.6M	125.6 g
Tris (free base)	121.14	1.2M	145.2 g
SDS (2%)	288.38	-	20.0 g
Sodium Meta-bisulfite	104.06	50mM	5.0 g
Ultrapure water (fill to)	-	-	1000 ml

* For non-reduced samples (especially antibodies), omit the Sodium Meta-bisulfite

Tris-Glycine-SDS Transfer Buffer, 10X, (Cat. # FB82500)

* pH should be between 8.4 and 8.6 at 25° C.

Ingredient	MW	Molarity	Qty/Liter
Tris (free base)	121.14	0.25M	30.3 g
SDS (2%)	288.38	-	20.0 g
Glycine	75.07	1.92M	144.1 g
Ultrapure water (fill to)	-	-	1000 ml

6.4 Recommended ClearPAGE™ Buffer Formulations, con't.

LDS Sample Buffer, 4x, (Cat. # FB31010)

*pH should be between 7.7 and 7.8 at 25° C.

Ingredient	MW	Molarity	Qty/Liter
Glycerol (40%)	-	-	400 g
Ficoll-400 (4%)	-	-	40 g
Triethanol amine, pH7.6	149.2	0.8M	120.0
6 N HCL	36.46	-	93.0 g
Lithium Dodecyl Sulfate (4%)	-	-	40 g
EDTA Di-Sodium	372.2	2mM	7.44 g
Brilliant Blue G250 (0.025%)	-	-	0.25g
Phenol Red	-	-	0.25 g
Ultrapure water (fill to)	-	-	1000 ml

Standard DNA/Native Running Buffer, 20X, (Cat. # GB61500)

*pH should be between 8.35 and 8.45 at 25° C.

Ingredient	MW	Molarity	Qty/Liter
Tricine (free acid)	179.17	0.8M	143.4 g
Tris (free base)	121.14	1.2M	145.2 g
Ultrapure water (fill to)	-	-	1000ml

DNA/Native Sample Buffer, 4x, (Cat. # GB33002)

* pH should be 7.6 at 25° C.

Ingredient	MW	Molarity	Qty/Liter
Glycerol (40%)	-	-	400 g
Ficoll-400 (4%)	-	-	40 g
Triethanol amine, pH7.6	149.2	0.8M	120.0
6 N HCL	36.46	-	93.0 g
EDTA Di-Sodium	372.2	2mM	7.44 g
Brilliant Blue G250 (0.025%)	-	-	0.25g
Phenol Red	-	-	0.25 g
Ultrapure water (fill to)	-	-	1000 ml

SECTION 7

Maintenance of Equipment

7.1 Care and Handling

The plastic components of the Quadra Electrophoresis Systems are fabricated from acrylic and polycarbonate. Electrodes and connectors are made from pure platinum, stainless steel, and chrome plated brass. As with any laboratory instrument, adequate care ensures consistent and reliable performance.

After each use, rinse buffer chamber, gel tray and combs with de-ionized water. Wipe dry with a soft cloth or paper towel, or allow to air dry. Whenever necessary, all components may be washed gently with water and a non-abrasive detergent, and rinsed and dried as above. Never use abrasive cleaners, glass cleaning sprays or scouring pads to clean the components, as these will damage the unit and components.

Additional precautions:

- Do not autoclave or dry-heat sterilize the apparatus or components.
- Do not expose the apparatus or components to phenol, acetone, benzene, halogenated hydrocarbon solvents or alcohols.
- Avoid prolonged exposure of the apparatus or components to UV light.
- Do NOT treat with diethylpyrocarbonate (DEPC)-treated water for extended periods at 37°C. A brief rinse with DEPC-water is sufficient after a thorough wash, followed by a quick rinse in 50% ethanol.

7.2 Maintenance

The following inspection and maintenance procedures will help maintain the safety and reliable performance of the Quadra Systems. Replacement parts can be ordered by calling 1-858-755-4959 or by contacting your local distributor.

- Banana plugs and power cords should be inspected regularly. If the banana plugs become loose or do not feel friction tight replace the plugs or power cords.
- Should power cord assemblies (connectors, wire or shrouds) show any signs of wear or damage (e.g. cracks, nicks, abrasions, or melted insulation), replace them immediately.
- The platinum wire is secured to the banana jack by compression between a stainless washer and the jack nut. The nut/washer interface should be tight and free of corrosion.

SECTION 8

Ordering Information

Quadra Ordering Information

COMBS: See page 6.

CAT. #	DESCRIPTION
QUADRA COMBO MINI-VERTICALS - 10CM(W) X 10CM(H)*	
QNX-700	Quadra System, CE, 4-place , Complete System
QNC-700	Quadra System, CE, 4-place, Cooled , Complete System
Accessories and Related Products for Quadra Systems	
EBX-BC-700	Accessory Blotting Cassette, includes sponge pad
EBX-SP-700	Accessory Sponge Pads, set of 4
QDC-P	Quick-Disconnect Coupling, self-closing
DCX-SP-108	Accessory shim plate for 8cm(w) Bio-Rad gels
DCX-P1010R	Gel Wrap® Glass Plate Set
DCX-E0710	Gel Wrap® Gasket, 0.75mm thick
DCX-E1010	Gel Wrap® Gasket, 1.0mm thick
DCX-E1510	Gel Wrap® Gasket, 1.5mm thick
DCX-S0710R	Gel Wrap® Spacer Set, 0.75mm thick
DCX-S1010R	Gel Wrap® Spacer Set, 1.0mm thick
DCX-S1510R	Gel Wrap® Spacer Set, 1.5mm thick
DCX-MC-4	Mini-Vertical Multi-Caster, 4-place, Unit Only
DCX-P1010	Multi-Caster Glass Plate Set
DCX-S0710	Multi-Caster Spacer Set, 0.75mm thick
DCX-S1010	Multi-Caster Spacer Set, 1.0mm thick
DCX-S1510	Multi-Caster Spacer Set, 1.5mm thick
EBU-204	Mini-Blotting System, CE, 4 place , fits 10 x 10cm gels
EBX-700	Mini-Blotting System, CE, 4 place , fits 10 x 10cm gels
EBU-4000	Semi-Dry Blotting System , transfer area of 20cm x 20cm
EPS-300 X	Power Supply, CE , 10-300V, 4-500mA
EPS-200 X	High Current Power Supply, CE , 5-200V, 10-2000mA

Quadra Multi-Caster Kit, 4-place - Cat. #'s

# Comb Wells	Spacer & Comb Thickness 0.75mm	Spacer & Comb Thickness 1.0mm	Spacer & Comb Thickness 1.5mm
10	DCX4K-0710	DCX4K-1010	DCX4K-1510
14	DCX4K-0714	DCX4K-1014	DCX4K-1514
16	DCX4K-0716	DCX4K-1016	DCX4K-1516

Quadra Cool Gel Wrap® Casting Bundle, 4-Place - Cat. #'s

# Comb Wells	Spacer & Comb Thickness 0.75mm	Spacer & Comb Thickness 1.0mm	Spacer & Comb Thickness 1.5mm
10	QN-C0710	QN-C1010	QN-C1510
14	QN-C0714	QN-C1014	QN-C1514
16	QN-C0716	QN-C1016	QN-C1516

Precast Gel Selection Chart & Ordering information

ClearPAGE SDS gels, packages of 10 - Cat. #'s					
Acrylamide %	Gel Size (cm) W x H	1 well Cat. #	2D well trough + 1 marker lane Cat. #	12 well Cat. #	17 well Cat. #
8%	10 x 8	-	-	BK00812-10	BK00817-10
	10 x10	-	-	FK00812-10	FK00817-10
10%	10 x 8	-	-	BK01012-10	BK01017-10
	10 x10	FK01001-10*	-	FK01012-10	FK01017-10
12%	10 x 8	-	-	BK01212-10	BK01217-10
	10 x10	FK01201-10*	-	FK01212-10	FK01217-10
16%	10 x 8	-	-	BK01612-10	BK01617-10
	10 x10	-	-	FK01612-10	FK01617-10
4 - 8%	10 x 8	-	-	BK40812-10	BK40817-10
	10 x10	-	-	FK40812-10	FK40817-10
4 - 12%	10 x 8	-	-	BK41212-10	BK41217-10
	10 x10	-	-	FK41212-10	FK41217-10
4 - 20%	10 x 8	-	-	BK42012-10	BK42017-10
	10 x10	FK42001-10*	FK42002-10*	FK42012-10	FK42017-10

CAT. #	DESCRIPTION
ClearPAGE™ BUFFERS, STAINS & MARKERS	

FB30001	Antioxidant, 1ml
FB30010	Antioxidant, 10ml
FB31010	LDS Sample Buffer, 4X, 10 ml
FB32001	DTT reducer, 10X, 1 ml
FB50053	SDS-non-reducing running buffer, 20X, 50 ml
FB50500	SDS-non-reducing running buffer, 20X, 0.5L
FB60053	Tricine SDS reducing running buffer, 20X, 50mls
FB60500	Tricine SDS reducing running buffer, 20X, 0.5L
FB14053	Turbo buffer SDS reducing running buffer, 20X, 50ml
FB14500	Turbo buffer SDS reducing running buffer, 20X, 0.5L
FB82500	Transfer buffer, with SDS, 10X, 0.5L
FB82600	Transfer buffer, without SDS, 10X, 0.5L
HG73010	Instant Blue stain, 1 L
HG73050	Instant Blue stain, trial size, 50ml
HM05160	2-color marker, 0.6ml
HM05103	2-color marker, trial size, 30µl

ClearPAGE™ BLOTTING AND DRYING	
FB82500	TGS transfer buffer, 10x, 500ml
HP19001	Blot sandwich NC 90 x 85mm, 0.2µl nitrocellulose/1mm blotting paper, 2 pack
HP19020	Blot sandwich NC 90 x 85mm, 0.2µl nitrocellulose/1mm blotting paper, 20 pack
HP29301	Blot sandwich PVDF 90 x 85mm, 0.2µl PVDF/1mm blotting paper, 2 pack
HP29320	Blot sandwich PVDF 90 x 85mm, 0.2µl PVDF/1mm blotting paper, 20 pack
HB04510	Gel drying solution, 1L
HB02520	Precut cellophane, 25cm x 25cm, 20 sheets
HB90001	Drying frame with 8 clamps with cellophane, 25cm x 25cm, 20 pieces of cellophane

ClearPAGE™ SAMPLER KITS	
FM90002	ClearPAGE Sampler Kit, Includes 2 each 10cm x 10cm SDS gels of customer's choice, trial buffers and instructions.
FM90005	ClearPAGE Super Sampler Kit, Includes 2 each 10cm x 10cm SDS gels of customer's choice, trial buffers, trial size 2-color marker, trial size Instant Blue stain and instructions.
FM90006	ClearPAGE Super Sampler Kit, Includes 2 each 10cm x 8cm SDS gels of customer's choice, trial buffers, trial size 2-color marker, trial size Instant Blue stain and instructions.

CONTACT INFORMATION

Online Ordering

www.cbsscientific.com

Telephone

Local and International: 858-755-4959

Toll Free: 800-243-4959

Sales E-mail Address

sales@cbssci.com

Technical Service E-mail Address

technicalservice@cbssci.com

Mailing Address

P.O. Box 856

Del Mar, CA 92014

Shipping Address

10805 Vista Sorrento Parkway, Ste 100

San Diego, CA 92121