

Selecting the best AbSelect™ Kit

Is your antibody purified?

Yes

Does the buffer contain BSA or Gelatin?

No

Does the buffer contain Tris, Glycine or Azide?

No

Is the antibody at least 1mg/ml?

Yes

No

AbSelect™ Antibody Concentration and Clean-up kit
[861-0010](#)

Yes

Yes

BSA only

AbSelect™ BSA Removal Kit
[820-0100](#)

Yes

Gelatin with or without BSA

AbPure™ Magnetic Purification System
Protein A Magnetic Beads
[265-0200](#)

AbSelect™ Antibody Purification Systems
(species optimised)

Rabbit
Protein A
[860-0010](#)

Human, sheep, goat, mouse, rat
Protein G
[890-0010](#)

Mouse
[830-0010](#)

No

What is the starting material?

Small volumes of Serum (less than 0.1ml) or Ascites Fluid (less than 0.5ml)

Large volumes of Serum or Ascites Fluid

Tissue culture supernatant (more than 10 ml)

AbSelect™ Serum Purification Systems
(species optimised)

Rabbit
Protein A
[863-0500](#)

Human, sheep, rat, goat, mouse
Protein G
[893-0030](#)

AbSelect™ TCS Purification Systems
(species optimised)

Rabbit, human, mouse, rat
Protein A
[862-0500](#)

Mouse
[832-0500](#)

Rat
[842-0500](#)

[Lightning-Link®](#), [Lightning-Link® Rapid](#) and [Thunder-Link® PLUS](#) Labeling kits